

HealthSky Embedding Cassettes

B3002-WH	75 pcs/sleeve	Lids Separate	3000 per Carton	White <input type="checkbox"/>
B3002-PK	75 pcs/sleeve	Lids Separate	3000 per Carton	Pink <input type="checkbox"/>
B3002-YL	75 pcs/sleeve	Lids Separate	3000 per Carton	Yellow <input type="checkbox"/>
B3002-RD	75 pcs/sleeve	Lids Separate	3000 per Carton	Red <input type="checkbox"/>
B3002-BR	75 pcs/sleeve	Lids Separate	3000 per Carton	Brown <input type="checkbox"/>

For more requirements, please contact HealthSky sales.

HealthSky Color Frosted Slides

B101000-WH	50 pcs/PK	Ground Edges 45° Corners	White <input type="checkbox"/>
B101000-PK	50 pcs/PK	Ground Edges 45° Corners	Pink <input type="checkbox"/>
B101000-GR	50 pcs/PK	Ground Edges 45° Corners	Green <input type="checkbox"/>
B101000-YL	50 pcs/PK	Ground Edges 45° Corners	Yellow <input type="checkbox"/>
B101000-BL	50 pcs/PK	Ground Edges 45° Corners	Blue <input type="checkbox"/>

For more requirements, please contact HealthSky sales.

SUPERIOR
PERFORMANCE
EVERY PRINT



EMARKER

Laser Printers

EMarker Laser Printers are a series of printing devices that deliver efficient and quality printing, enable workflow optimization.

Hangzhou HealthSky Biotechnology Co., Ltd.

Tel/Fax: +86 571 89774886

E-mail: overseas@healthskybio.com Web: www.healthskybio.com

Add.: Floor 3, Building 2-1, Liangzhu International Life Science and Technology Town, Yuhang District, Hangzhou, Zhejiang, China.

HEALTHSKY



HEALTHSKY
Histology Solutions

About HealthSky

Founded in 2004, HealthSky has been dedicated to the research, development, manufacturing, sales, and service in the field of pathology for two decades.

HealthSky is committed to providing unique and advanced solutions for cervical cancer screening and histologic diagnosis.

Located in Hangzhou, HealthSky spans an area of 12,000 square meters and is certified with ISO 13485:2016. The facility includes an R&D center, a standard manufacturing plant, and an accredited pathology diagnostic laboratory. The company employs over 300 people, including more than 50 R&D researchers, with 60% holding Ph.D. or Master's degrees.

With a deep understanding of pathology, a professional R&D team, and state-of-the-art facilities, HealthSky is revolutionizing the process of pathologic diagnosis to make it more standardized, automated, intelligent, and precise.

12 Invention Patents

35 Utility Patents

4 Software Copyrights



Pathology, the bridge between science and medicine. HealthSky will keep striving for excellence, working for a healthier and brighter future for humans.



HealthSky Template Editor

HealthSky Template Editor is used to design cassette or slide templates using texts, serial numbers, icons, and barcodes.

It allows users to create the design from the start or by editing the template by freely adding, deleting or moving the elements as well as changing their properties.



User-friendly Interface



Freely Template Design



Easy LIS Integration



Flexible Connectivity



HealthSky Print Asst.

EMarker PrintAsst. manages the running of the printer, as well as operation and data processing.

It is used to select templates, edit the printing data, sequence, and perform daily print tasks.

EBMarker-160

Multi Hopper Laser Cassette Printer



EMarker-21A-C600

Multi Hopper Ultraviolet Laser Cassette Printer

EMarker-21A-C600 utilizes laser printing technology to create crisp, permanent, and high-resolution markings at high speed.

It is a six-hopper cassette printer with a large capacity and supports hopper selection printing to meet the cross-selection of different tissue cassette types.



Model No.	EMarker-21A-C600	EBMarker-160
Printing Technology	Ultraviolet Laser	Infrared Laser
Compatible Cassette	45° angle marking surface, with or without lid.	POM material with laser absorber, 45° angle marking surface, with or without lid.
Print Speed	3~4 Second/ Cassette	
Resolution	2500 dpi	
Load Capacity	6 Hoppers	
Print Content	Alphanumeric, 1D and 2D barcodes, Graphics, Logos	
Auxiliary Functions	Repetitive task reminder, Running fault alarm	
Operation & Display	PC	
Software Compatibility	Windows11, Windows10, Windows 7, x64/x86; Windows XP	
Data Interface	Wire	
Integration	LIS/HIS	
Supply Voltage	100-240 VAC, 50/60 Hz	
Dimension (L × W × H)	320×420×510 (mm)	
Weight	24 kg	22 kg



Superior Print Quality



High Print Speed



Freely Hopper Selection

EMarker-21A-S100

Ultraviolet Laser Slide Printer

EMarker-21A-S100 utilizes laser printing technology that can leave permanent marks on the writing surface of the microscope slides.



EMarker-11A-S100

Laser Slide Printer

EMarker-21A-C100

Ultraviolet Laser Cassette Printer

EMarker-21A-C100 is a laser printer that enables efficient, dependable, and error-free printing of cassettes.



EMarker-11A-C100

Laser Cassette Printer



Economical



Practical



Flexible



Compact Size



Scan to Print



Convenient Slide Collection

Model No.	EMarker-21A-S100	EMarker-11A-S100
Printing Technology	Ultraviolet Laser	Infrared Laser
Compatible Slide	Color frosted slides	Color frosted slides with laser absorber
Print Speed	3~4 Second/ Slide	
Resolution	2500 dpi	
Load Capacity	100 Slides	
Print Content	Alphanumeric, 1D and 2D barcodes, Graphics, Logos	
Auxiliary Functions	Dust filter, Scan input, Automated slide collector(optional), Repetitive task reminder, Running fault alarm	
Operation & Display	PC	
Software Compatibility	Windows11, Windows10, Windows 7, x64/x86; Windows XP	
Data Interface	Wire	
Integration	LIS/HIS	
Supply Voltage	100-240 VAC, 50/60 Hz	
Dimension (L × W × H)	245×410×360 (mm)	220×420×380 (mm)
Weight	18.5 kg	15 kg

Model No.	EMarker-21A-C100	EMarker-11A-C100
Printing Technology	Ultraviolet Laser	Infrared Laser
Compatible Cassette	45° angle marking surface, with or without lid	POM material with laser absorber, 45° angle marking surface, with or without lid.
Print Speed	3~4 Second/ Cassette	
Resolution	2500 dpi	
Load Capacity	1 Hopper	
Print Content	Alphanumeric, 1D and 2D barcodes, Graphics, Logos	
Auxiliary Functions	Repetitive task reminder, Running fault alarm	
Operation & Display	PC	
Software Compatibility	Windows11, Windows10, Windows 7, x64/x86; Windows XP	
Data Interface	Wire	
Integration	LIS/HIS	
Supply Voltage	100-240 VAC, 50/60 Hz	
Dimension (L × W × H)	280×410×360 (mm)	220×420×380 (mm)
Weight	18.5 kg	15 kg